

I. Preparing and titering baculovirus working stocks for protein production.

Because we can only send you a small sample of a plaque-purified, first passage (PP1P1) virus stock, you will need to produce and titer a working virus stock before producing your recombinant protein of interest. You should grow the virus by infecting Sf9 or Sf21 cells cultivated in TN-MFH medium containing 10% FCS and 0.1% pluronic F68 (“complete” TN-MFH). Please NOTE that baculovirus stock produced in serum-free media are not stable. For more detailed protocols, see (Jarvis, 2011; Murhammer, 2007; O'Reilly, Miller, and Luckow, 1992), among others.

A) Amplify the PP1P1 stock:

- 1) Aseptically seed 0.5×10^6 Sf9 cells/ml of complete TN-MFH medium into a shake flask.
- 2) Incubate overnight in a shaking incubator at 28°C, 125 rpm.
- 3) Aseptically add ~0.1 mL of the recombinant baculovirus PP1P1 stock.
- 4) Incubate in a shaking incubator at 28°C until you see clear signs of infection:
 - a) Look for swollen cells, swollen nuclei.
 - b) This usually takes ~3-5 days.
- 5) Aseptically harvest the medium and pellet cells and cell debris:
 - a) Ten min at top speed (~1000 g) in a tabletop centrifuge.
- 6) Aseptically remove the supernatant, transfer to a fresh tube or bottle:
 - a) This is the PP1P2 virus stock.
 - b) Wrap with aluminum foil to protect from light and store at 4°C (Jarvis and Garcia, 1994).
 - (i) These stocks will slowly lose titer during storage (~1 log/yr; Jarvis and Garcia, 1994).
 - (ii) Can freeze aliquots for longer-term storage/archiving.

B) Titer the working virus stock:

- 1) Aseptically seed 0.75×10^6 Sf9 cells/well into one or more 6-well plates:
 - a) Cells will stick best if pelleted, resuspended in serum-free TN-MFH.
 - b) Use 3 ml per well for seeding.
 - c) Allow cells to adhere for ~1 h at 28°C.
- 2) Aseptically prepare a 10-fold dilution series (10^{-1} to 10^{-6}) with PP1P2 stock:
 - a) Use complete TNM-FH as the diluent.
 - b) Can typically plate (10^{-4} , 10^{-5} , 10^{-6}) and get a useful range of plaques to count.
- 3) Drain medium from cells.
- 4) Replace with 2 mL of diluted virus in each well:
 - a) Plate each dilution in duplicate.
 - b) Incubate for ~1 h at 28°C.
- 5) Prepare agarose overlay:
 - a) Melt 1.5 mL of 2% agarose/well in a microwave oven.
 - b) Equilibrate in a 60°C water bath.
 - c) Warm 1.5 mL of complete 2X Grace's medium/well in an 30°C water bath.
- 6) Overlay the infected cell monolayers:
 - a) Aseptically mix equal volumes of agarose and complete 2X Grace's medium (3.0 mL/well).
 - b) Mix thoroughly.
 - c) Aseptically remove the virus from each well in the plaque assay.
 - d) Overlay with 3 mL per well.
 - e) Allow agarose to harden for 10 min.
- 7) Seal in a plastic baggie and incubate upside down for 7-10 days at 28°C.
- 8) Remove from baggie, count plaques on upside down plate under a dissecting microscope.
 - a) Titer = avg # plaques X reciprocal of the dilution (pfu/mL).

II. Protein production and purification.

A) Infect insect cells to produce the protein of interest.

- 1) Aseptically seed 0.5×10^6 Sf9 cells/ml in your favorite insect cell medium into a shake flask:
 - a) Best to use Sf9 cells adapted to a serum-free medium to purify a secreted protein.
- 2) Incubate in the shaking incubator overnight at 28°C.
- 3) Perform a viable cell count to determine cell density (should ~double overnight).
- 4) Aseptically add the volume of virus needed to achieve a multiplicity of infection of 2-5 pfu/cell.
- 5) Incubate in a shaker-incubator at 28°C for 2-3 days:
 - a) Typically, 48-60 h.
 - b) But, best to perform preliminary experiments to determine the time of infection that yields the highest levels of intact product for each recombinant protein.
- 6) Harvest the medium and pellet cells and cell debris:
 - a) Ten min at top speed (~1000 xg) in a tabletop centrifuge.
- 7) Remove the supernatant, place in a fresh tube or bottle as the cell-free fraction.
 - a) For secreted proteins, process the supernatant as the source of recombinant protein.
 - b) For intracellular proteins, process the cellular pellet as the source of recombinant protein.
- 8) For detailed information on immobilized metal affinity chromatography, see:
 - a) (Crowe et al., 1994).
 - b) (Petty, 2001).
 - c) Vendors instructions.
 - d) Many papers on HIS-tagged proteins expressed in baculovirus-infected insect cells.

B) Metal affinity-purify the protein of interest (secreted proteins).

- 1) Dialyze the supernatant against dialysis buffer (50 mM NaH_2PO_4 , 500 mM NaCl; pH 8.0) for 12-18 h with a buffer change after 6 h.
- 2) Wash Probond (Invitrogen) or your favorite metal affinity resin with dialysis buffer.
- 3) Incubate the dialyzed supernatant with the washed resin at RT for 3 h with gentle agitation.
- 4) Transfer into an empty column barrel and allow the resin to settle into a bed.
- 5) Wash with 5 column volumes of dialysis buffer containing low concentration of imidazole:
 - a) We use 20-50 mM imidazole for washes.
- 6) Elute with 1 column volume of dialysis buffer containing high concentration of imidazole:
 - a) We use 250 mM imidazole.
- 7) Assess purification by analysis of starting material, unbound, wash, and elution fractions.

C) Metal affinity-purify the protein of interest (intracellular proteins).

- 1) Wash the infected cells once or twice with ice cold TBS (50 mM Tris, 100 mM NaCl; pH 8.0).
- 2) Lyse the cells in one of two ways:
 - a) Resuspend pellet in a minimal volume of ice-cold dialysis buffer; freeze-thaw twice.
 - b) Resuspend pellet in a minimal volume of TBS + 1% NP40; incubate for 10 min on ice.
- 3) Shear DNA by titration with an 18-gauge needle until the solution is no longer viscous.
- 4) Centrifuge for 15 min at ~10,000 xg at 4°C.
- 5) Harvest the supernatant and use for metal affinity chromatography as described above (D2).

III. References.

- Crowe, J., Dobeli, H., Gentz, R., Hochuli, E., Stuber, D., and Henco, K. (1994). 6xHis-Ni-NTA chromatography as a superior technique in recombinant protein expression/purification. *Methods Mol Biol* 31, 371-87.
- Jarvis, D. L. (2011). Recombinant protein expression in baculovirus-infected insect cells. In "Methods Navigator Protocols (online)". Elsevier, Amsterdam.
- Jarvis, D. L., and Garcia, A., Jr. (1994). Long-term stability of baculoviruses stored under various conditions. *Biotechniques* 16(3), 508-513.
- Murhammer, D. W. (2007). "Baculovirus and insect cell expression protocols." 2nd ed. Methods in Molecular Biology, 388 Humana Press.
- O'Reilly, D. R., Miller, L. K., and Luckow, V. A. (1992). "Baculovirus expression vectors." W.H. Freeman and Company, New York.
- Petty, K. J. (2001). Metal-chelate affinity chromatography. *Curr Protoc Mol Biol* Chapter 10, Unit 10 11B.