



PCR Quantification with Melt Curve Report PCR Base Line Subtracted Curve Fit Data Graphs Only


Current Date: **11-Oct-06 12:54 PM**
Data generated on: **12-Sep-05 at 10:09 AM.**

Optical data file name: **Val 91205.odm**
Plate Setup file used: **5fold4samples.psm**
Protocol file used: **Alison78DC.tmo**

Sample volume: **20.00 ul**
Hot Start? **No**
Well factor collection: **Persistent (Experimental Plate Failed)**

Comments


PCR Amp/Cycle Graph for SYBR-490



Data Analysis Parameters

Calculated threshold has been replaced by the user selected threshold **50.0**.
Per-well baseline cycles have been determined automatically.
Data analysis window is set at **95.00%** of a cycle, centered at **end** of the cycle.
Weighted Mean digital filtering has been applied. Global filtering is **off**.

Melt Curve Graph for SYBR-490



Melt Curve Analysis Parameters

Weighted Mean digital filtering has been applied. Global filtering is **off**.
Threshold for automatic peak detection is set at **1.00**.